



Problem solver

FitLife provides comfort for hard-to-fit patients.

The FitLife total face mask is a solution for patients who have difficulty wearing conventional masks or who simply prefer a clear line of sight.

Facial hair, dentures, nasal bridge breakdown, facial irregularities, and claustrophobia are issues that can greatly impact the ability to establish effective therapy and comfort. FitLife offers a solution for men and women facing these challenges. It quickly establishes a robust seal around the less pressure-sensitive perimeter of the face, while maintaining an unrestricted view. Easy-to-use snap clips and advanced headgear make application and removal straightforward. The large size of this “total solutions mask” fits most patients. The new FitLife is the latest way that we help you fit 100% of your patients, 100% of the time.

PHILIPS

RESPIRONICS

FitLife headgear is available in small and large sizes.

EZ Peel tabs allow for quick, easy adjustment and removal of headgear.

Snap clips simplify re-application and eliminate the need to refit after removal.

Over-molded cushion is integrated into the mask frame, reducing the number of parts.

Perimeter seal with soft silicone cushion fits quickly and evenly for maximum comfort.

Large surface area equalizes pressure and flow inside the mask, minimizing eye irritation.

Entrainment valve triggers access to room air if pressure drops below 3 cm H₂O.



Accessory swivel allows for freer tubing movement.



Built-in radial diffusion leak ports deflect exhaled air away from the bed partner and provide quiet operation.

In-line pressure pick-off port allows for pressure measurement or O₂ bleed in.

Item	Part number	HCPCS code
FitLife mask		
Small mask with small headgear	1060801	A7030, A7035
Large mask with large headgear	1060802	A7030, A7035
Parts and accessories	Part number	HCPCS code
FitLife headgear, small	1060807	A7035
FitLife headgear, large	1060808	A7035
FitLife headgear snap clips (2 per pack)	1060913	N/A
Accessory swivel	7041	N/A

Respironics, FitLife, and EZ Peel are trademarks of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Philips Healthcare is part of Royal Philips Electronics

Europe, Middle East, Africa:
+33 1 47 52 30 00
Asia Pacific:
+852 3194 2280

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668
+1 724 387 4000
+1 800 345 6443